

CYCV14-CHANGYU pTOPO-TA Kit

Contents	20 rxns (CYCV1401)	80 rxns (CYCV1402)
pTOPO-T Vector (30ng/μl)	20 μl	80 μl
1000bp Control (30ng/μl)	5 μl	5 μl
10 × Enhancer	20 μl	80 μl

Note: All reagents, when stored at -20 °C, are stable for 12 months.

Description

The Zero Background TA Topoisomerase Cloning Kit is designed for fast cloning of DNA fragments up to 10 kb generated by *Taq* DNA polymerase which introduces single 3'-A overhang to the end of the DNA product. DNA fragments generated by restriction digestion or obtained by mechanical shearing can also be effectively cloned after end polishing using *Taq* DNA polymerase. It utilizes DNA strand transfer activity of *Vaccinia* virus topoisomerase I. *Vaccinia* virus DNA topoisomerase I forms a 3'-phosphoryl intermediate with the plasmid vector containing cleavage recognition motif of 5'CCCTT ↓. Covalently bound topoisomerase I then transfers the incised vector DNA strand to the DNA fragment to be cloned with free 5'-OH termini. This transferring reaction is rapid and reproducible. The cloning vector pTOPO-TA included in this kit is a high-copy-number plasmid engineered to tolerate mildly toxic genes. Regions flanking the cloning site of pTOPO-TA vectors have multiple common restriction sites for release of the cloned fragment by single or double restriction digestion.

❖ Procedure

1. Set up the topoisomerase cloning reaction by mixing the reagents in the order shown.

DNA Fragment or 1μl 1000bp control	0.5-8μl
pTOPO-T Vector	1μl
10 × Enhancer	1μl
ddH ₂ O	Xμl
Final Volume	10μl

DNA fragments MUST NOT be 5'-phosphorylated. TA-ended DNA fragment is preferred. Highest cloning efficiency was observed from PCR fragments generated by *Taq* DNA polymerase or fragments tailored by *Taq* polymerase. 10-150 ng of

DNA fragments ranging from 100 bp to 5000 bp have been tested to give satisfactory results. Excessive amounts of DNA inserts (e.g., >300 ng) will significantly reduce cloning efficiency.


Fragment size (bp)	Optimal amount (ng)
100-1000	10-40
1000-2000	40-80
2000-5000	80-150

- Mix the reaction gently and incubate for 5 minutes at room temperature between 20-30°C.

Recent R&D data indicate incubation at this step produces no detectable benefits. In any case, do not let the incubation go beyond 5 minutes. Extended incubation for larger inserts up to 5 kb is unnecessary and may introduce background.

- Add 5µl of the cloning reaction into 50 µl chemically competent E. coli and mix gently. Do not mix by pipetting up and down.
- Incubate on ice for 2-30 minutes.
- Follow the competent cell manufacturer's instructions to complete the transformation.

❖ **Map of pTOPO-T:**




❖ **Sequencing Primer:**

M13F: TGTAACACGACGGCCAGT

M13R: CAGGAAACAGCTATGACC

❖ **MCS information:**



This product is furnished for **LABORATORY RESEARCH USE ONLY**.
 Not for diagnostic or therapeutic use.